







Actual

size

Press one button and the concentration of Soy Milk is displayed within seconds!

within seconds! Just run it under the faucet.

## Features

- Requires only a few drops of sample
- Wait 3 seconds and it's done!
- Compact: Fits right in your Palm!
- Only weighs 100g (less than half a pound!)
- Automatic Temperature Compensation (ATC)
- Works indoors and outdoors with "ELI" function
- Hang it from your neck or keep it in your pocket!
- Can be washed under the faucet! (IP65)
- Convenient storage case
- Battery power indicator
- Designed with a strap holder option

## **Specifications**

Measurement	Concentration of Soy Milk
Measurement range	0.0 to 20.0%
Resolution	0.1%
Measurement accuracy	±0.2%
Measurement temperature	10 to 60°C (ATC)
Ambient temperature	10 to 40°C
Measurement time	3 seconds
Power supply	AAA alkaline battery × 2
Battery life	Approx. 11,000 measurements (Using new alkaline batteries)
International Protection Class	IP65 Water Resistant
Dimensions and weight	55(W)×31(D)×109(H)mm, 100g

到 (0~20

Pocket

ZERO

PAL-27S

HACCF GMP All ATAGO refractometers are designed and manufactured in Japan.

## ATAGO CO., LTD.

Headquarters: The Front Tower Shiba Koen, 23rd Floor 2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan TEL : 81-3-3431-1943 FAX : 81-3-3431-1945 overseas@atago.net http : //www.atago.net/

- CATAGO U.S.A., Inc. TEL : 1-425-637-2107 customerservice@atago-usa.com
- CATAGO INDIA Instruments Pvt.Ltd. TEL:91-22-2833-8038 / 8076 customerservice@atago-india.com
- CATAGO (THAILAND) Co., Ltd. TEL : 662-982-8718-9 customerservice@atago-thailand.com
- CATAGO BRASIL Ltda. TEL : 55 16 3916-6000 customerservice@atago-brasil.com
- ATAGO ITALIA s.r.I. TEL : 39 2 36557267 customerservice @atago-italia.com CATAGO CHINA Guangzhou Co., Ltd. TEL : 86 20-38106065 info@atago-china.com



\* Specifications and appearance are subject to change without notice.

GLP

ATAGO products comply with HACCP, GMP, and GLP system standards.